

Instrument and service trolleys



Instrument trolley 609

Frame with chromed steel tube. Two detachable stainless steel plates (40 x 60 cm). 75 mm castors make the trolley easy to move.

The trolley can be placed under a 80 cm high worktop.



Instrument trolley 611, includes a drawer

Frame with chromed steel tube. Two detachable stainless steel plates (40 x 60 cm). 75 mm castors make the trolley easy to move.

Underneath the upper plate there's an easy-to-clean, plastic accessory box that can be opened on both sides. The trolley can be placed under a 80 cm high worktop.



Ward trolley 262

Frame with chromed steel tube and a push handle. Two solid stainless steel plates (40 x 60 cm). Big 100 mm castors make the trolley easy to move.

Includes an accessory box of easy-to-clean plastic. The box can be opened on both sides.



Service trolley 266 and Courier trolley 265

Frame with chromed steel tube and a push handle. Baskets (2pcs) are epoxy-coated steel wire. Big 100 mm castors.

The upper basket in Courier trolley 265 is 6-sectional.



Service trolley 285

Frame is stainless steel. Two solid stainless steel plates (55 x 80 cm) with edges to prevent items from falling . Push handle. Big 125 mm castors (2 lockable) of corrosion resistant material.



Samples trolley 255

Frame with chromed steel tube. Push handle. Detachable plastic boxes with edges to prevent items from falling. Castors 75 mm.



Linen trolley 270

Frame with chromed steel tube. Four twin baskets of epoxy-coated steel wire, designed to fit both mechanically and manually folded linen. Castors 100 mm.

	Instrument trolley 609	Instrument trolley 611	Ward trolley 262	Service trolley 266	Courier trolley 265	Service trolley 285	Samples trolley 255	Linen trolley 270
Length	600 mm	600 mm		942 mm	942 mm	800 mm	620 mm	990 mm
Height	750 mm	750 mm	780mm	750 mm	750 mm	900 mm	750 mm	1400 mm
Width								505 mm
Castors	75	75	100	100	100	125	75	100

Soiled linen trolleys



SOILED LINEN TROLLEY 290

Trolley for two sacks. Base supports and sack fastening mechanism. Frame with chromed steel tube. Castors Ø 75 mm. Designed for 100-150 liter fabric or plastic bags.

Accessories:

Lid



SOILED LINEN TROLLEY 291

Trolley for one sack. Base support and sack fastening mechanism. Lid. Push handle. Castors Ø 75 mm. Designed for 100-150 liter fabric or plastic bags.

Accessories:

Garbage bag holder for a 30 liter garbage bag



SOILED LINEN TROLLEY 292

Trolley for one sack. Base support, sack fastening mechanism and foot pedal. Lid. Push handle. Castors Ø 125 mm. Designed for 100-150 liter fabric or plastic bags.