

ANTIMICROBIAL UPHOLSTERIES DISINFECTANT RECOMMENDATION

Manufacturer	Disinfectant	Silvery
Antibac AS	Antibac	XX
B. Braun Melsungen AG	Meliseptol foam pureHexaquart	XX
Berner Pro	Perasafe	XX
Berner Pro	Virkon	XX
BODE Chemie GmbH	Kohrsolin	XX
Desomed Dr. Trippen GmbH	Desomed rapid	XX
Dr. Schumacher GmbH	Cleanisept Wipes	XX
Dürr Dental AG	FD360	XX
Dürr Dental AG	FD360	XX
Ecolab	Brial	XX
Ecolab	Klercide 60/40	XX
Ecolab	Incidin Oxywipe	XX
KiiltoClean Oy	Easydes	XX
KiiltoClean Oy	Klorilli	XX
Laboratoires Anios	Surfanios	XX
Laboratoires Anios	Surfa ´ safe	XX
Lysoform Dr. Hans Rosemann GmbH	Aldasan	x
Schülke & Mayr GmbH	Antifect plus	x
Schülke & Mayr GmbH	Buraton rapid	-
Schülke & Mayr GmbH	Microzid	XX
Wet Wipe A/S	Wet wipe	XX

xx = recommended

x = recommended to use diluted, 0.50% dilution

- = not recommended